

PARTICLE SIZE ANALYSIS SET: PIPETTE METHOD

Used for determining the particle size distribution of very fine soils.

Basic items to perform the test are the 22-T0062/1 Andreasen pipette which is hold by the 22-T0062/2 Pipette stand for raising or lowering precisely the pipette to the required level without disturbing the suspension under test. The 22-T0062/3 sedimentation cylinder is submerged in the 22-T0058/A Constant temperature water bath.

22-T0058/A, /AZ

- Dimensions: 600x300x380 mm
- Weight approx.: 10.5 kg

22-T0062/2A

- Weight approx.: 10 kg



PARTICLE SIZE ANALYSIS SET: PIPETTE METHOD

STANDARDS

- BS 1377:2

ORDERING CODE

22-T0058/A

Constant temperature water bath. 230V/50Hz/Ph

22-T0058/AZ

Constant temperature water bath. 110V/60Hz/Ph

22-T0062/1

Andreasen pipette

22-T0062/2A

Stand for precisely positioning (up-down) the Andersen pipette to the required height without disturbing the test.

